HOPE Car Bed


Unique travel solutions for children with special health care needs requiring a larger car bed

Weight Range 4.5 - 35 lbs.



Merritt Manufacturing Inc.
Listening to Therapists; Working for Families ™

Large Car Bed Features



© 2024, Merritt Manufacturing, Inc.

Features

Five Restraint Options

Medical personnel to determine position and type of restraint to be used



Supine Position
Restraint Bag with
adjustable
cummerbund



Prone Position
Restraint Bag with
adjustable cummerbund



Side Facing Position
(Right Side only)
Restraint Bag with
adjustable cummerbund



Supine Position
3-point harness with adjustable cummerbund



Prone Position
3-point harness with adjustable cummerbund



Child Restraint Bag

with adjustable cummerbund

two sizes available



Restraint Bag features dual entry zippers

Bag securely attaches to bed



Cummerbund available in two sizes -- adjusts for best position on child



Prone position in restraint bag must have medical approval



Child Restraint Bag

- Standard car bed equipped with small restraint bag /small cummerbund and 3-point harness/large cummerbund to cover full range of weights (4.5 – 35 pounds.)
- Large restraint bag available for children over 10-pounds to 35 pounds.
- Small and large side facing restraint bags available.



Small Restraint Bag 4.5 to 10 pounds

Small and Large Side Facing Restraint Bags

Large Restraint Bag 10 to 35 pounds



3-Point Harness - Supine Position

Easy to adjust 3-point Harness for children 10-35 pounds in supine position



Cummerbund to position child and protect from buckle and chest clip

Buckle can be positioned to desired location



3-Point Harness

3-Point Harness - Prone Position

Prone position must have medical approval




Cummerbund adjusts for best position on child in prone position.

Protects child from buckle and chest clip while in prone position.

- 3-point Harness option for prone position for a child 10 to 35 lbs.
- Soft padded adjustable cummerbund is required for positioning and comfort



Adjustable buckle with vinyl pad provides best fit for child





Spica Cast Options



Modified Small Restraint Bag Child 4.5 - 10 pounds un-casted

3 pt Harness Child over 10 pounds un-casted



© 2024, Merritt Manufacturing, Inc.

Installation in the vehicle

For Vehicles with rear lap/shoulder belts

To be installed in rear seat on the right hand side behind the front seat passenger position. The right hand and center position lap-shoulder belts are to be used....

Or with right hand lap-shoulder belt and center position lap belt.

Can also be installed using the right hand lap belt and center lap belt.

Left seating position lap-shoulder belt is for the parent or

medical attendant.

Note: Some smaller vehicles may not allow room for the parent or medical attendant.

Child's head must be positioned toward the center of the vehicle to protect the child from side impact.





Installation in the vehicle

Leveling Straps

Includes two car bed leveling straps if front seat frames are available

This will provide four attachment points and will be the most secure installation



Foam seat wedge available for leveling of the car bed on vehicle seats with more severe slopes or if leveling straps are not used







Hope Car Bed Features

- Largest Car Bed available with 5 Restraint options
- Designed for single patient, hospital loaner programs and medical institutional use
- Child's weight range: 4.5 35 lbs.
- Child's Height: Up to 29" (longer if legs are permitted to bend)
- Bed size: Inside dimensions 14.25" wide, 29.5" long
- Car bed weight: 22.2 lbs. Includes handles for easy transport outside the vehicle
- 6 year plastic shell life -- all parts replaceable
- Easy to install securely in vehicle



Hope Car Bed Questions



Contact:

Merritt Manufacturing Inc. 580 Tracy Rd. New Whiteland, IN 46184-9699

> Phone: 317-409-0148 Fax: 317-893-2567

www.merrittcarseat.com

E-mail: info@merrittcarseat.com

Made in USA

Designed & manufactured in the United States by Merritt Manufacturing, Inc.