

# CHURCHILL SPECIFICATIONS

Minimum weight: **44** lbs.  
Maximum weight: 175 lbs.

Minimum height: **44** inches  
Maximum height: 72 inches

Shoulder Width Maximum: **Unlimited**  
Thigh Width Maximum: **Unlimited**

Seat depths:  
– 15”  
– 18”

Super low profile seat simplifies transfer of occupant into vehicle.

First belt positioning booster for people with special health care needs weighing between **44-175** lbs.

Ask us how we can customize your Churchill!

As always, we accept returns and charge NO restocking fees!

See DEMO VIDEOS:  
[www.youtube.com/user/MerrittManufacturing](http://www.youtube.com/user/MerrittManufacturing)



# merritt

listening to therapists | working for families™

Designed and manufactured in the United States by Merritt Manufacturing, Inc., The Churchill™ conforms to all applicable Federal motor vehicle safety standards.

©2024 Merritt Manufacturing, Inc. All Rights Reserved

## MERRITT MANUFACTURING, INC.

580 Tracy Rd.  
New Whiteland, IN 46184-9699  
Phone: 317-409-0148  
Fax: 317-893-2567  
Email: [info@merrittcarseat.com](mailto:info@merrittcarseat.com)

[www.MerrittCarSeat.com](http://www.MerrittCarSeat.com)  
Check us out on social media!



**THE CHURCHILL™**  
A belt-positioning booster seat specially designed for individuals with special health care needs.



# merritt

listening to therapists | working for families™

# FEATURES

\*DENOTES STANDARD FEATURE

## POSITIONING HARNESS\*



For complete torso support, the vest is fully adjustable and easy to put on. The vest is supported by two vertical straps. Available in three sizes.

## BUCKLE GUARD



Prevents the child from opening their car seat buckle.

## CHEST CLIP GUARD



Prevents the child from pushing car seat chest clip down or opening it.

## "EZ-UP" HEAD REST SYSTEM



Developed for occupants with poor head control, the EZ-Up™ system enables the occupant to maintain a natural head position. Use with the EZ-Up Cap™ to restrict forward motion.

## "STAY-PUT" POMMEL



This pommel (abductor) stays firmly in place while in use. The pommel pivots away for easy transfer.

## HIP INFLECTION WEDGE



Provides anti-thrust by raising the occupant's knees by 2.5 inches.