



**P.O. Number:**

**Date:**

## 2024 Chamberlain Order Form

**The Chamberlain is available in four fabric color options. To select a color, simply fill in the space at the end of the part number with the color option abbreviation from the following list.**

- Sapphire Blue – S
- Rapunzel Purple - RP
- Total Camo – TC
- Beige Naugahyde – N

**See Important Statement Below About Installing the Chamberlain in vehicles manufactured prior to 2003\*\***

**Billing Address**

Contact Name:
Organization:
Street Address:
City, State:
ZIP Code:
Phone: <span style="float: right;">Fax:</span>
E-Mail:

**Shipping Address (If Different)**

Contact Name:
Organization:
Street Address:
City, State:
ZIP Code:
Phone: <span style="float: right;">Fax:</span>
E-Mail:

Qty	Part Number	Description	Unit Price	Line Total
	4000CHSS-___	<b>The Chamberlain Standard with Small Positioning Vest</b> – an adaptive pediatric and adult positioning vehicle system (Weight range: 81 to 225 lbs). Small Positioning Vest - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches.		
	4000CHSM-___	<b>The Chamberlain Standard with Medium Positioning Vest</b> – an adaptive pediatric and adult positioning vehicle system (Weight range: 81 to 225 lbs). Medium Positioning Vest - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches.		
	4000CHSL-___	<b>The Chamberlain Standard with Large Positioning Vest</b> – an adaptive pediatric and adult positioning vehicle system (Weight range: 81 to 225 lbs). Large Positioning Vest - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches.		
	<b>Please</b>	<b>One vest is included in the price of</b>		

	<b>Note*:</b>	<b>the Chamberlain. Pricing for additional vests is as follows.</b>		
	2000CVS-___	<b>Small Positioning Vest - for chest circumference between 25 and 32 inches and belly circumference less than 44 inches. *Price charged only when ordering additional vests.</b>		
	2000CVM-___	<b>Medium Positioning Vest - for chest circumference between 30 and 37 inches and belly circumference less than 44 inches. *Price charged only when ordering additional vests.</b>		
	2000CVL-___	<b>Large Positioning Vest - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. *Price charged only when ordering additional vests.</b>		
	2000VB	<b>Add Side Release Buckles to one of the vests above. (Side Release buckles can help stop an occupant from escaping the vest).</b>		
		<b>Optional Accessories</b>		
	2000EUHR-___	EZ-Up Head Rest – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below:		
	<b>See note to the right regarding the five cap color options.</b>	<b>Note: There are five options for cap color. The Khaki option is the coolest option for hot temperatures. Simply fill in the space at the end of the part number below with “K” for Khaki or one of the color abbreviations (to match your seat cover) at the top of this order form.</b>		
	1000EUCS-___	EZ-Up Cap – Size Small - (Head Circumference: 17.5 – 19.5 inches)		
	1000EUCM-___	EZ-Up Cap – Size Medium - (Head Circumference: 19.5 – 21.5 inches)		
	1000EUCL-___	EZ-Up Cap – Size Large - (Head Circumference: 21 – 23.25 inches)		
	2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – <b>See Below**</b>		
			<b>Sub Total***</b>	
<b>Indiana Residents Only; Not Applicable to other States</b>			<b>Sales Tax (7%)</b>	
			<b>TOTAL</b>	

\* The “Chamberlain Standard” includes the following features:

- Torso Support System
- Anti-Slouching Leg Strap Assembly
- Positioning Vest (Small, Medium, or Large)
- Tether Extension Assembly

**\*\* The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Chamberlain must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 – 2002 models also have LATCH (See your vehicle owner’s manual to verify if LATCH is present).**

**For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle’s floor board.**

**Please call Customer Service at 317-409-0148 with any questions regarding the installation of the Chamberlain.**

**\*\*\* Shipping charges – The cost of shipping “Best Way” will be added to Purchase Orders.**

**Please place orders by one of the following methods:**

**E-Mail: [order@merrittcarseat.com](mailto:order@merrittcarseat.com)**

**Phone: 317-409-0148**

**Fax: 317-893-2567**

**Thanks for Your Business!**

**Merritt Manufacturing, Inc.  
[www.merrittcarseat.com](http://www.merrittcarseat.com)**

**580 Tracy Rd,  
New Whiteland, IN 46184**