



P.O. Number:

Date:

2024 Churchill with Positioning Vest Order Form

The Churchill-Pediatric Positioning
Vehicle Restraint System and Accessories

See Important Statement on Last Page About Installing the Churchill****

The Churchill is available in four fabric color options. To select a color, simply fill in the space at the end of the part number with the color option abbreviation from the following list.

- Sapphire Blue – S
- Rapunzel Purple - RP
- Total Camo – TC
- Beige Naugahyde - N

Billing Address

| |
|--|
| Contact Name: |
| Organization: |
| Street Address: |
| City, State: |
| ZIP Code: |
| Phone: Fax: |
| E-Mail: |

Shipping Address (If Different)

| |
|--|
| Contact Name: |
| Organization: |
| Street Address: |
| City, State: |
| ZIP Code: |
| Phone: Fax: |
| E-Mail: |

| Qty | Part Number | Description | Unit Price | Line Total |
|-----|-------------|--|------------|------------|
| | 2000CS15-__ | The Churchill Standard* with Vest and 15 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). One Positioning Vest is included in price. Please choose vest size below. | | |
| | 2000CS18-__ | The Churchill Standard* with Vest and 18 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). One Positioning Vest is included in price. Please choose vest size below. | | |

| | | | |
|---|---|--|--|
| 2000CVS-___ | Small Positioning Vest - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. **Price charged only when ordering vest separately. | | |
| 2000CVM-___ | Medium Positioning Vest - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. **Price charged only when ordering vest separately. | | |
| 2000CVL-___ | Large Positioning Vest - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. **Price charged only when ordering vest separately. | | |
| 2000CVC-___ | Custom Vest – Please Provide Following Measurements: Chest Circumference: ____; Waist Circumference: ____; Height from Hip to Shoulder: ____ | | |
| 2000VB | Add Side Release Buckles to one of the vests above. (Side release buckles can help stop an occupant from escaping the vest. For the proficient escape artist, rather than the vest with buckles, we suggest selecting the Churchill with Positioning Harness and adding the Buckle guard and the Chest clip guard. See "Churchill with Positioning Harness" in the product menu). | | |
| | Optional Accessories | | |
| 2000EUHR-___ | EZ-Up Head Rest – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest. | | |
| See note to the right regarding the five cap color options. | Note: There are five options for cap color. The Khaki option is the coolest option for hot temperatures. Simply fill in the space at the end of the part number below with "K" for Khaki or one of the color abbreviations (to match your seat cover) at the top of this form. | | |
| 1000EUCS-___ | EZ-Up Cap – Size Small - (Head Circumference: 17.5 – 19.5 inches) | | |
| 1000EUCM-___ | EZ-Up Cap – Size Medium - (Head Circumference: 19.5 – 21.5 inches) | | |
| 1000EUCL-___ | EZ-Up Cap – Size Large - (Head Circumference: 21 – 23.25 inches) | | |
| 2000FA-___ | Foam Abductor | | |
| 2000HICP-___ | Hip Inflection Comfort Pommel | | |
| 2000HIW-___ | Hip Inflection Wedge | | |
| 2000IC-15-___ | Quick-Change Incontinent Cover for 15 Inch Seat Depth Base | | |

| | | | | |
|---|---------------|--|-----------------------|--|
| | 2000IC-18-___ | Quick-Change Incontinent Cover for 18 Inch Seat Depth Base | | |
| | 2000ADD | Adductor Strap | | |
| | 2000AB | Attachment Bracket Kit (for Vehicles without LATCH) – See Important Note Below**** | | |
| | | | | |
| | | | Sub Total** | |
| Indiana Residents Only; Not Applicable to other States | | | Sales Tax (7%) | |
| | | | TOTAL | |

* The "Churchill with Positioning Vest" includes the following features:

- Base
- Lower Anchor Attachment Straps
- Anti-Slouching Leg Strap Assembly
- Positioning Vest
- Foam Padding
- Fabric Cover
- Tether Extension Assembly
- Belt Positioning Guides

** One Vest is included in price of "Churchill with Positioning Vest."
Retail Price of additional vests is \$375.00 each.

**** The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 – 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board. Please call Customer Service at 317-409-0148 with any questions regarding the installation of the Churchill.

Please place orders by one of the following methods:

E-Mail: order@merrittcarseat.com
Phone: 317-409-0148
Fax: 317-893-2567

Thanks for Your Business!

** Shipping charges – The cost of shipping
"Best Way" will be added to Purchase Orders.

Merritt Manufacturing, Inc.,
580 Tracy Rd,
New Whiteland, IN 46184

www.merrittcarseat.com