



## 2024 Churchill with Positioning Harness - Retail Price List

The Churchill-Pediatric Positioning  
Vehicle Restraint System and Accessories

**See Important Statement Below About Installing the Churchill\*\*\*\***

**The Churchill is available in four fabric color options. To select a color, simply replace the asterisk in the part numbers in the table with the color option abbreviation from the following list.**

- Sapphire Blue – S
- Rapunzel Purple - RP
- Total Camo – TC
- Beige Naugahyde - N

Part Number	Description	Retail Price
2000CH15-*	<b>The Churchill with Positioning Harness*** and 15 Inch Seat Depth - <u>CAN BE USED WITH OUR ESCAPE PROOF PRODUCTS</u></b> – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). <b>One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.</b>	\$995.00
2000CH18-*	<b>The Churchill with Positioning Harness*** and 18 Inch Seat Depth - <u>CAN BE USED WITH OUR ESCAPE PROOF PRODUCTS</u></b> – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). <b>One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.</b>	\$995.00
	<b>Optional Accessories</b>	
2000BG	Buckle Guard - (To stop occupant from opening buckle)	\$150.00
2000CCG-*	Chest Clip Guard - (To stop occupant from opening Chest Clip or pushing it down)	\$150.00
<b>Important Note</b>	<b>One of the following positioning vests can be added for use with the position harness. This would also allow the possibility of using only the vest; or only the harness.</b>	

2000CVS-*	Small Positioning Vest - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches.	\$375.00
2000CVM-*	Medium Positioning Vest - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches.	\$375.00
2000CVL-*	Large Positioning Vest - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches.	\$375.00
2000CVC-*	Custom Vest – Please Provide Following Measurements: Chest Circumference: ____; Waist Circumference: ____; Height from Hip to Shoulder: ____	\$375.00
2000VB	Add Side Release Buckles to one of the vests above.	\$99.00
2000EUHR-*	EZ-Up Head Rest – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest.	\$195.00
See note to the right regarding the five cap color options.	<b>Note: There are five options for cap color. The Khaki option is the coolest option for hot temperatures. Simply fill in the space at the end of the part number below with "K" for Khaki or one of the color abbreviations (to match your seat cover) at the top of this price list.</b>	
1000EUCS-*	EZ-Up Cap – Size Small - (Head Circumference: 17.5 – 19.5 inches)	\$125.00
1000EUCM-*	EZ-Up Cap – Size Medium (Sapphire Blue) - (Head Circumference: 19.5 – 21.5 inches)	\$125.00
1000EUCL-*	EZ-Up Cap – Size Large (Sapphire Blue) - (Head Circumference: 21 – 23.25 inches)	\$125.00
2000FA-*	Foam Abductor	\$125.00
2000HICP-*	Hip Inflection Comfort Pommel	\$175.00
2000HIW-*	Hip Inflection Wedge	\$105.00
2000IC15-*	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base	\$145.00
2000IC18-*	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base	\$145.00
2000ADD	Adductor Strap	\$75.00
2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – <b>See Important Note Below****</b>	\$95.00

**\*\*\*The "Churchill with Harness" includes the following features:**

- Base
- Lower Anchor Attachment Straps
- Positioning Harness
- Foam Padding
- Fabric Cover
- Tether Extension Assembly
- Belt Positioning Guides

**\*\*\*\* The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 – 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).**

**For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.**

**Please call Customer Service at 317-409-0148 with any questions regarding the installation of the Churchill.**

Shipping charges – The cost of shipping "Best Way" will be added to Purchase Order

**Please place orders by one of the following methods:**

**E-Mail: [order@merrittcarseat.com](mailto:order@merrittcarseat.com)**

**Phone: 317-409-0148**

**Fax: 317-893-2567**

**Thanks for Your Business!**

**Merritt Manufacturing, Inc.  
[www.merrittcarseat.com](http://www.merrittcarseat.com)**

**580 Tracy Rd,  
New Whiteland, IN 46184**